

IDEXX xChek*: Creating a User-Defined Assay for the IDEXX Influenza A Virus Antibody Test—Avian Species

Creating a user-defined assay is an easy process of two tasks—editing an xChek* file to allow you to create an assay, and then creating the assay within xChek.


To edit the xChek.ini file:

1. From the task bar, choose **Start> Run**, and type “xchek.ini” into the text box. Click the **OK** button.
2. After the xChek.ini file opens, locate the line labeled “AllowNewAssays.”
3. If the value for that line is “0,” change it to “1.” If the value is already “1,” do not do anything.
4. Locate the section labeled [BasicGraphs], scroll to the end of the section, and then enter “InfAavian=17.”
5. From the menu bar, choose **File> Save**, and then choose **File> Exit**.


To create the assay:

1. From the task bar, choose **Start> Programs> xChek> xChek** to open the xChek program.
2. From the menu bar, choose **Database> Assays> All** to open the Assay dialog box.
3. Click the **New** button, and enter “Influenza A Virus Antibody Test Kit Avian Species” into the Name field, “InfAavian” into the Code field, and “unknown” into the Kit Lot field. Today’s date (mm/dd/yyyy) is automatically entered into the Expiration field.

Note: When you begin testing, update the information in the Kit Lot and Expiration fields with the actual kit lot number and expiration date.



4. Click the **Basic** button, and select or enter the following information for the respective fields:
Case Type: Population Samp Filter: 650
Template: Poultry Vert Ref. Filter: 0
Species: Chicken Dilution: 10
Technology: ELISA Wells: 1



5. Click **OK** to save the settings.

- Click the **Calculations** button, and select or enter the following information:

Formulas:

"S" Part of Ratio: Sample1

"P" Part of Ratio: None

"N" Part of Ratio: Negative

Blocking Factor: None

First Calculation:

Variable: S/N

Positive Cutoff: 0.5

Suspect Cutoff: 0.5

Comparison: <

Second Calculation:

Variable: None

Positive Cutoff: 0

Suspect Cutoff: 0

Comparison: None

The image shows a software dialog box titled "New Assay". It is divided into several sections. The "Formulas" section contains four dropdown menus: "S" Part of Ratio (set to "Sample1"), "P" Part of Ratio (set to "None"), "N" Part of Ratio (set to "Negative"), and Blocking Factor (set to "None"). To the right of these are three buttons: "OK", "Cancel", and "Reset". Below this are two columns of settings. The "First Calculation" column has four dropdowns: Variable (set to "S/N"), Positive Cutoff (set to "0.5"), Suspect Cutoff (set to "0.5"), and Comparison (set to "<"). The "Second Calculation" column has four dropdowns: Variable (set to "None"), Positive Cutoff (set to "0"), Suspect Cutoff (set to "0"), and Comparison (set to "None").

- Click **OK** to save these settings, and then click the **OK** button to save the assay to the database. You can now test for influenza A antibodies using the xChek Assay Management System.

Note: xChek does not evaluate controls or results for validity when you use a user-defined assay. You must evaluate the results from each assay in accordance with good laboratory practices. To evaluate your assay validity, refer to the insert provided with the test kit.

Setting Up Additional Display Variables

Only OD values are displayed in the reporting options for this user-defined assay. To set up additional display variables, follow the directions below.

- From the xChek menu bar, choose **Reports> Analyze Cases** to open the Filter Criteria for Analyze dialog box.
- Enter the desired search criteria and click **OK**. The Select Cases for Analyze dialog box appears.

3. Select the desired cases and click **OK**. The Analyze Report dialog box appears.

Analyze Report

Case Options

- Show Block Reports
- Show Statistics
- Show Controls
- Show NHC Wells
- Show Original ODs
- Show Comments
- Show Kit Lot/Expiration
- Single Column

Graph Options

- Show Graphs
- Show Statistics
- Show Key
- Color Graphs
- Stats on Top

Mean

- Arithmetic
- Geometric

Style

- 2D Histogram
- 3D Histogram

OK

Cancel

Variables

Sort

Baseline

Footers

4. Click the **Variables** button to open the Display Variables dialog box.

Display Variables

- Well Type
- Optical Density
- Tube Number
- Sex
- Species
- Breed 1
- Breed 2
- Animal Type
- Animal ID1
- Animal ID2
- Calendar Age
- Breeding Age
- Quality
- Mean Sample
- Sample
- Pos
- Neg
- Blocking
- S/P Ratio
- S/N Ratio
- B/A Ratio
- S/NHC Ratio
- Titer
- Titer Log2
- Titer Group
- Result

OK

Cancel

5. Select the desired options and click **OK**. Click **OK** again to save the settings.

For more information, call IDEXX Technical Services at 00800 727 43399, contact your local area manager, or visit us on the Web at idexx.com.

Test With Confidence™

IDEXX

*xChek and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. in the United States and/or other countries.

© 2011 IDEXX Laboratories, Inc. All rights reserved. • 9540-00

One IDEXX Drive
Westbrook, Maine 04092 USA